DR/F System Retrofit

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Spec #	Spec Item
C-1	1024 × 1024 pixel resolution
C-2	12-bit digitization depth
C-3	Data digitization at camera's end for improved Signal to Noise Ratio
C-4	Computerized TV camera iris control
C-5	Comprehensive Fluoroscopy Acquisition Package
C-6	Fluoroscopy Loop Acquisition
C-7	Simultaneous acquisition and processing
C-8	PCI- based Workstation
C-9	Windows- Based Mouse-Driven User Interface
C-10	Automatic windowing, storage and filming
C-11	Radiography acquisition rate of up to 7.5 frames per second
C-12	Real-time disk
C-13	Anatomically programmed protocols
C-14	User-selectable gamma corrections
C-15	Off-line image subtraction
C-16	Full Compliance with DICOM 3.0 Standard

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- Sensor : CCD Frame Transfer
- Resolution : 1024 * 1024 , Progressive and interlace scan
- Dynamic range : 1000 :1
- Digitization depth : 12bit
- Antiblooming : greater than x100
- Frame rates : Up to 30 fps in progressive mode
- Up to 60 fps in interlaced mode
- Sampling rate : 40MHz
- Lens with computer controlled iris
- Automatic gain control (AGC)

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(2) Automatic gain(Brightness) control

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(1) Patient review

User-selectable images (studies/series) may be transfer-red to other applications for review / processing.

(2) Window / level

Continuously variable window limits of the displayed image. The operator may program six preset windows.

(3) Window polarity inversion

Positive/negative polarity inversion

(4) Image inversion

Enable rotation of the image by 90,180, or 270 deg. Horizontal and Vertical reverses are available as well.

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대항목	소항목	특징
Video Inputs	4 Video Inputs	0.5 V _{pp} ~ 2V _{pp} , 75 Ω termination or open-loop, 12bit gain, 12bit BW level
	RGB Input	3-pass
	Sync Source	Composite sync on any analog channel
	Bandwidth	150MHz @ 1.0 Vpp input
	Horizontal Frequency	Up to 105KHz
	External Clock	1.0 Vpp ~ 5.0 Vpp, up to 150MHz
Video Format	RS-170, RS-343, RS-330, CCIR, non-standard video signals	
Imaging Quality	S/N Ratio	59dB @ 50 MHz
	Linearity	Better than 1%
	Gain and Offset Stability	1%
	Pixel Jitter	±0.5ns
	Synchronization Time	Less than 250 ms
Programmable Characteristics	Horizontal Resolution	Up to 4,096 total pixel clocks per line
	Horizontal Resolution	Up to 4,096 total lines
	Horizontal Capture Delay	12 bit register with resolution of 1 pixel
	Vertical Capture Delay	12 bit register with resolution of 1 line
Acquisition Characteristics	A/D Conversion	10 bits per pixel: no missing codes
	Capture Method	Single pass up to 75 MHz; dual pass up to 150 MHz
Memory Characteristics	Control Registers	Memory mapped: PCI Plug 'n' Play
	Storage Memory 4 Megapixels	Data Transfer Speed PCI bus master: bursts up to 132 MB/sec: 45 MB/sec sustained
Software	Supports up to 4 boards per system Windows 98, Windows 95, Windows NT 4.0 Auto-SYNC automatic configuration software and example programs are included with each HI*DEF Accura. The IDEA Software Developer's Kit can be purchased separately. Drivers for Image-Pro Plus, Optimas, and TWAIN are available.	

Interactive zoom, with magnification factor of up to X4.

(6) Magnifying glass

Zoom of X2 in a user-selectable Region of interest.

(7) Alpha numeric-text annotation

User-selectable text can be superimposed anywhere in the image. Text may be moved / removed within the image.

(8) Graphic aids

Including are arrows, lines and Rols. Distance measurement is also available.

(9) General image filtration (up to 7x7 kernel).

72 preset edge filters (II size dependent) as well as pressersmooth filters are available.

(10) Cine loops with user-controllable display rates

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(11) Electronic shutter

Data outside the image's active area is blackened, thus increasing image conspicuity.

Determination of the active area is performed either automatically or under user's control.

(12) Gamma collections

Different gamma correction curves (user-selectable) are available for the in-room monitor, The workstation monitor and the laser imager.

(13) Multiformat display

Formats of 2,4,6,9,12 and 20 images within a single frame.

(14) Image editing Copy, Duplicate, Cut,Paste, Swap, Delete

(15) Automatic of all post processing operations is also available.

Images are transferred to the printer in the background, thus Enabling the performance of further acquisition/ processing tasks.

(16) Optimized imaging parameters

Several imaging parameters are automatically set by the system.

According to the acquisition technique/anatomical program.

Selected by the operator.

(17) Options (Angiography package)

① Real-time subtraction (DSA)

- The mask image is subtracted on-line from all acquired images.
- Real-time subtracted images are displayed during sequence acquisition.

② Road-mapping

- Two-phases protocol for accurate placement of catheters during Fluoroscopic angiography procedures.
- During the first phase, a mask image is generated on-line from the fluoroscopic Images using a minimum intensity projection algorithm.
- During phase II of the protocol, the mask is subtracted on-line from all incoming Averaged images.

③ Pixel shifting

- The operator can modify the relative position between the mask and contrast images.
- Shifting resolution are 1/8,1/4,1/2 and 1 pixel.

- Pixel shifting may be performed Manually or automatically.
- ④ DSA with maximum/minimum pixel opacification.
 - An optimum DSA image is obtained using maximum or minimum Intensity projection algorithm.
- ⑤ Group Mask / Group Subtraction
 - Optimized mask and contrast images are obtained using averaging techniques.
 - Images are integrated with a 16-bit accuracy.
- ⑥ Land marking
 - Mask with variable weighing factor is created and subtracted from all selected Images this enabling the display of the background anatomy.

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- (1) Dicom 3.0 storage service class
- (2) Dicom 3.0 printer service class
- (3) Dicom 3.0 modality Work-list service class