

# *In vitro* study of fracture strength of provisional crown materials

#### Isil Karaokutan<sup>1\*</sup>, Gulsum Sayin<sup>2</sup>, Ozlem Kara<sup>3</sup>

<sup>1</sup>Department of Prosthodontics, Faculty of Dentistry, Selcuk University, Konya, Turkey <sup>2</sup>Department of Prosthodontics, Faculty of Dentistry, Medipol University, Istanbul, Turkey <sup>3</sup>Department of Prosthodontics, Faculty of Dentistry, Bezmialem Vakif University, Istanbul, Turkey

**PURPOSE.** The purpose of this report was to evaluate the effect of the fabrication method and material type on the fracture strength of provisional crowns. **MATERIALS AND METHODS.** A master model with one crown (maxillary left second premolar) was manufactured from Cr-Co alloy. The master model was scanned, and the data set was transferred to a CAD/CAM unit (Yenamak D50, Yenadent Ltd, Istanbul, Turkey) for the Cercon Base group. For the other groups, temporary crowns were produced by direct fabrication methods (Imident, Temdent, Structur Premium, Takilon, Systemp c&b II, and Acrytemp). The specimens were subjected to water storage at 37°C for 24 hours, and then they were thermocycled (TC, 5000×, 5-55°C) (n=10). The maximum force at fracture (Fmax) was measured in a universal test machine at 1 mm/min. Data was analyzed by non-parametric statistics ( $\alpha$ =.05). **RESULTS.** Fmax values varied between 711.09-1392.1 N. In the PMMA groups, Takilon showed the lowest values (711.09 N), and Cercon Base showed the highest values (959.59 N). In the composite groups, Structur Premium showed the highest values (1392.1 N), and Acrytemp showed the lowest values (910.05 N). The composite groups showed significantly higher values than the PMMA groups (*P*=.01). **CONCLUSION.** Composite-based materials showed significantly higher fracture strengths than PMMA-based materials. The CAD-CAM technique offers more advantages than the direct technique. [J Adv Prosthodont 2015;7:27-31]

**KEY WORDS:** Temporary crown; Provisional restoration; Composite; Methyl methacrylate; Fracture strength, CAD-CAM

# **INTRODUCTION**

Provisional restorations are an important part of prosthetic therapy procedures with fixed prostheses (i.e., crowns and bridges).<sup>1</sup> Provisional restorations serve important roles during tooth preparation and until fitting, luting the final fixed restoration.<sup>1-3</sup> These include pulpal tissue protection against physical, chemical, and thermal injuries; mainte-

© 2015 The Korean Academy of Prosthodontics This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (http://creativecommons. org/licenses/by-nc/3.0) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited. nance of positional stability and occlusal function; and provision of the prepared teeth with strength, retention, and aesthetics, which are essential to clinical success. Polymethyl methacrylate (PMMA) resins and composite-based resins (CBR) are the most common materials used to fabricate provisional fixed dental prosthesis (FDP).<sup>2,4,5</sup> Their chemical natures differ; methacrylate resins use liquid/powder and are hand-mixed, and composite-based resins use paste/ paste and are usually auto-mixed. The polymerization reaction of methacrylate resins initiates chemically (self-curing), while composite-based materials are available as both selfcuring and dual-curing systems.

Fractures are a common cause of failure of provisional restorations. Although restorations should be designed to avoid failure, fractures can still occur. This may cause the patient discomfort and economic loss. Thus, the mechanical strength properties of provisional materials are important and should be considered to ensure the clinical success of provisional restorations.<sup>6</sup> Incorrect occlusion, bruxism, undercontoured pontics, and trauma are potential reasons

Corresponding author:

Isil Karaokutan

Department of Prosthodontics, Faculty of Dentistry, Selcuk University, Semsi Tebrizi Mah. Ankara Cad. No:6 42030 Karatay-Konya, Turkey Tel. 90 332 223 1202: e-mail, dt. isilsentoregil@gmail.com Received May 20, 2014 / Last Revision September 22, 2014 / Accepted September 30, 2014

for restoration fractures during usage.<sup>3,7</sup> Fractures can occur even during normal masticatory functions, especially when the patient has long-span bridges.<sup>8,9</sup>

Chemically polymerized materials available for provisional restorations using either PMMA or CBR have unique physical properties that are dependent upon the composition of the chemical monomer. Different monomers have different effects such as exothermic reactions, polymerization shrinkage, marginal fit, periodontal responses, color stability, and fracture strength.<sup>10,11</sup> Fracture strength is related to the mechanical properties of the provisional restorative materials.<sup>9,11,12</sup> Previous studies on the resins used for provisional restorations compared the mechanical properties of PMMA and CBR; however, the results were controversial.<sup>12-14</sup>

Chairside-fabricated temporary restorations are associated with some shortcomings regarding mechanical strength, surface texture, and fit<sup>15-17</sup>; certain mixing procedures and overfilling of the impression might lead to voids that compromise their mechanical strength.<sup>16</sup> Furthermore, after producing these restorations, the flexural strength showed very low values.<sup>15</sup>

CAD/CAM technologies, which are used to fabricate temporary restorations, may solve some of these issues; resin-based blanks cured under optimum conditions exhibited increased mechanical strength and prevented porosities within the restorations.<sup>18</sup> In addition, CAD/CAM-fabricated temporary restorations reportedly reduced the chairside time and produced superior results.<sup>19</sup>

This study investigated the influence of different fabrication methods and materials on the fracture strength of different provisional crowns.

# MATERIALS AND METHODS

The fracture resistance of different materials and fabrication methods was evaluated in a laboratory assembly (25°C, 50% rel. humidity) on a Cr-Co alloy master model with a crown FDP. Table 1 shows the materials tested in this study. All specimens were prepared according to the manufacturers' recommendations. A non-carious human upper left second premolar (no. 25) was used in this study. The tooth was embedded in a chemically cured acrylic resin block to 2 mm below the cementoenamel junction, and it was duplicated and cast using a Cr-Co alloy (Brealloy C+ B 270, Bredent GmbH, Senden, Germany). Shoulder preparation (convergence angle 6°) was used for the full crown preparation. The prepared tooth was also duplicated and cast using a Cr-Co alloy.

Six types of provisional restorative materials were directly fabricated using the over-impression technique. A vinyl polysiloxane impression material (Imprint<sup>TM</sup> 3 VPS Impression Material, 3M ESPE, MN, USA) was mixed, and it was placed in the casted master model. After setting, the impression was removed and controlled for any damage to the contour. Afterwards, the temporary crown was mixed according to the manufacturer's recommendations and placed into the over impression. After setting according to the manufacturer's recommendation, the over impression was cut into two pieces, and the temporary crown material was carefully removed. This procedure was used for all six groups (n=10).

The master model was scanned, and the data set was transferred to a CAD/CAM unit (Yenamak D50, Yenadent Ltd, Istanbul, Turkey) for the Cercon Base group. The FDPs were milled and carefully removed from the blocks.

The specimens were controlled for the existence of air bubbles. Inaccurate specimens were discarded and replaced. Ten crowns per experimental group were subjected to deionized distilled water at 37°C for 24 hours, and then they were thermocycled for 1 week (TC, 5000×, 5-55°C; dwell time 30 seconds, transfer time 2 seconds).

The Cr-Co-alloy-prepared tooth model was placed in a universal testing machine (LS 500; Lloyd Instruments, West Sussex, UK) for fracture testing. The specimens were placed on a model. Specimens were subjected to a compressive load at a 90-degree angle to the center of the specimen until fracture at a crosshead speed of 1 mm/min Fmax was recorded. Fracture patterns were evaluated optically.

The influence of materials and fabrication methods was analyzed by the Kruskal-Wallis test (P<.05). Post hoc com-

#### Table 1. Materials used in this study

Product	Manufacturer	MR*	Main components of the monomer mixture	
Imident	Imicryl Dis Malzemeleri/Konya, Turkey	1:2.4	Polymethyl methacrylate	
Structur Premium	VOCO GmbH/Cuxhaven, Germany	1:1	Bisacrylic composite resin	
Systemp c&b ll	Ivoclar Vivadent AG/Schaan, Liechtenstein	4:1	Polyurethane polymethacrylate	
Acrytemp	Zhermack spA/Via Bovazecchino, Italy	4:1	Bisacrylic composite resin	
Takilon BBF	WP GmbH/Barmstedt, Germany	2.1/10	Polymethyl methacrylate	
Temdent Classic	Schütz Dental GmbH/Rosbach, Germany	NA	Polymethyl methacrylate	
Cercon Base	DuguDent GmbH/Hanau, Germany	NA	Highly cross-linked methyl methacrylate	

\* Mixing ratio dimethacrylates-base:catalyst [by volume]; mixing ratio monomethacrylate-liquid:powder [volume:mass].

parisons were calculated by Mann-Whitney tests (P<.05). All statistical analyses were performed using SPSS for Windows (12.0, SPSS Inc., Chicago, IL, USA).

## RESULTS

Both the materials and the fabrication techniques had a significant effect on the fracture strength (P<.05). Figure 1 shows the results of the descriptive analyses. As seen by the mean values of the CBR groups, Structur Premium showed the highest values (1392.1 ± 344.11 N), and Systemp c&b II showed the second highest (1009.0 ± 84.50 N). Acrytemp showed the lowest values of the CBR groups (910.05 ± 77.09). The PMMA groups Takilon (711.09 ± 179.18 N) and Temdent (745.23 ± 94.75 N) showed the lowest values. Imident showed the highest values of directly fabricated PMMA groups (843.71 ± 83.46). The only CAD-CAM PMMA group, Cercon Base, showed similar results (1106 ± 134.65 N) to the CBR Systemp c&b II group.

Table 2 shows the significance of the results of the Mann-Whitney U test. As shown, there is no significant difference between the CBR group Systemp c&b II and the

CAD-CAM PMMA Cercon Base groups (P=.07), nor between the PMMA group Imident and the CBR group Acrytemp (P=.123). The other PMMA groups, Takilon and Temdent, also showed no significant difference (P=.290). The CBR group Structur Premium and the CAD-CAM PMMA group Cercon Base showed some significance (P=.049). The PMMA group Imident showed significant differences from the Temdent, Systemp c&b II, Takilon, Structur Premium, and Cercon Base groups (P<.05). The CBR group Acrytemp showed significant differences from the Temdent, Systemp c&b II, Takilon, Structur Premium, and Cercon Base groups (P < .05). The PMMA group Temdent showed significant differences from the Systemp c&b II, Structur Premium, and Cercon Base groups (P<.05). The CBR group Systemp c&b II showed significant differences from the Takilon and Structur Premium groups (P<.05). The PMMA group Takilon showed significant differences from the Structur Premium and Cercon base groups (P < .05). The CAD-CAM PMMA group Cercon Base showed highly significant differences from other PMMA groups (P < .05).



Fig. 1. Descriptive analyses of materials.

**Table 2.** Mann-Whitney U test significance values (Significant at *P*<.05)

	Imident	Acrytemp	Temdent	Systemp c&b II	Takilon	Structur Premium	Cercon Base
Acrytemp	.123		<.01	.016	.003	<.01	<.01
Temdent	.034	<.01		<.01	.29	<.01	<.01
Systemp c&b II	<.01	.016	<.01		<.01	<.01	.07
Takilon	.19	<.01	.29	<.01		<.01	<.01
Structur Premium	<.01	<.01	<.01	<.01	<.01		.049
Cercon Base	<.01	<.01	<.01	.07	<.01	.049	

## DISCUSSION

In the present study, the fracture resistance of seven temporary crown materials was evaluated. While this *in vitro* study may not reflect the oral conditions, strength values could be a useful predictor of clinical performance and helpful for comparing provisional materials tested in a controlled situation.

A temporary luting cement was intentionally excluded to omit an additional influencing variable. It was assumed that the luting cement would have increased the fracture strength; this subject should be investigated in further studies.

Restorations become weaker when they are exposed to compatible temperature changes (5-55°C) similar to those found in an intraoral environment. In this study, all specimens were placed in distilled water at 37°C for 24 hours, and then were long-term thermocycled (5000×, 5-55°C, dwell time 30 seconds, transfer time 2 seconds).<sup>20,21</sup> This TC is equal to 6 months of clinical use. Some studies have investigated the effect of different storage times on temporary materials.<sup>12,22,23</sup> Balkenhol et al.<sup>15</sup> found a positive correlation between storage time and mechanical properties. Koumjian and Nimmo9 found similar results and further discovered that dry storage showed higher transverse strength values than wet storage for all materials. Conversely, some previous investigations have shown relatively large increases in the mechanical properties of certain bis-acryl and PMMA interim resin materials between 1 hour and 24 hours storage times.15,24,25,26

Different chemical compositions can be responsible for differences in fracture strength. Conventional methyl methacrylate-type resins are monofunctional; they have a low molecular weight and are linear molecules that exhibit decreased strength and rigidity.<sup>27</sup> Indeed, two of the materials that displayed the lowest fracture strengths were methacrylate resins (Temdent and Takilon).

Bis-acryl resin composite materials are difunctional, and thus they are capable of cross-linking with another monomer chain. This cross-linkage provides strength and durability to the material. Two of the materials with the highest fracture strengths were bis-acryl resins (Voco Structur and Systemp c&b II). The CAD-CAM group Cercon Base showed the highest fracture strength among PMMA groups and was higher than the Acrytemp and Systemp c&b II composite groups. The manufacturer of the Cercon Base PMMA stated that the material included highly cross-linked PMMA and was cured under idealized conditions. Thus, the CAD-CAM PMMA material is a more convenient temporary material than the other PMMA groups made by direct techniques.

The higher mechanical strength of acrylic-based temporary crowns compared to traditional monomethacrylates is in concurrence with the literature. Nejatidanesh *et al.*<sup>28</sup> found that bis-acryl provisional materials showed higher flexural strength than methacrylate resins. Lang *et al.*<sup>29</sup> compared two PMMA and four composite temporary materials in an artificial oral environment and found that the highest strength values were accompanied by low fracture rate in the composite-based group.

In this study, all the fractures were seen on the force applied surfaces and all the failures were recorded as catastrophic. Occlusal forces are usually measured in the incisors and molars in a stomatological system. The average values are 250 N in the incisor zone and 350 N in the molar zone, but the values are much higher in patients with bruxism.<sup>30</sup> Some studies have shown that these values can rise to 720-815 N.<sup>30,31</sup> Two of our PMMA groups (Temdent and Takilon) fractured within these values. For patients with bruxism, clinicians should choose the temporary crown and bridge material with care.

## **CONCLUSION**

According to the findings of this study, composite-based materials offer more advantages than PMMA-based materials in regard to fracture strength. Thus, they should be preferred as a material for provisional restorations. PMMA-based CAD/CAM fabricated provisional crowns show higher fracture strength then directly fabricated crowns. Computer aided design and manufacturing might increase the strength of provisional restorations. Further studies will be supported by a greater variety of temporary CAD/CAM materials.

### REFERENCES

- Shillingburg HT, Sumiya H, Lowell DW, Jacobi R, Brackett SE. Fundamentals of fixed prosthodontics. 3. Auflage. Chicago: Quintessence; 1997.
- Burns DR, Beck DA, Nelson SK; Committee on Research in Fixed Prosthodontics of the Academy of Fixed Prosthodontics. A review of selected dental literature on contemporary provisional fixed prosthodontic treatment: report of the Committee on Research in Fixed Prosthodontics of the Academy of Fixed Prosthodontics. J Prosthet Dent 2003;90: 474-97.
- Gough M. A review of temporary crowns and bridges. Dent Update 1994;21:203-7.
- Anusavice KJ. Phillips' science of dental materials. 10<sup>th</sup> ed. Philadelphia; WB Saunders; 1996; p. 285-431.
- Gegauff AG, Holloway JA, Rosenstiel SF, Land MF, Fujimoto J. Contemporary fixed prosthodontics. 3<sup>rd</sup> ed. Mosby, St. Louis; 2001; p. 380-416.
- Kim SH, Watts DC. In vitro study of edge-strength of provisional polymer-based crown and fixed partial denture materials. Dent Mater 2007;23:1570-3.
- Burke FJ, Murray MC, Shortall AC. Trends in indirect dentistry: 6. Provisional restorations, more than just a temporary. Dent Update 2005;32:443-4, 447-8, 450-2.
- Rosentritt M, Behr M, Lang R, Handel G. Flexural properties of prosthetic provisional polymers. Eur J Prosthodont Restor Dent 2004;12:75-9.
- 9. Koumjian JH, Nimmo A. Evaluation of fracture resistance

of resins used for provisional restorations. J Prosthet Dent 1990;64:654-7.

- Young HM, Smith CT, Morton D. Comparative in vitro evaluation of two provisional restorative materials. J Prosthet Dent 2001;85:129-32.
- 11. Diaz-Arnold AM, Dunne JT, Jones AH. Microhardness of provisional fixed prosthodontic materials. J Prosthet Dent 1999;82:525-8.
- Ireland MF, Dixon DL, Breeding LC, Ramp MH. In vitro mechanical property comparison of four resins used for fabrication of provisional fixed restorations. J Prosthet Dent 1998;80:158-62.
- 13. Osman YI, Owen CP. Flexural strength of provisional restorative materials. J Prosthet Dent 1993;70:94-6.
- Wang RL, Moore BK, Goodacre CJ, Swartz ML, Andres CJ. A comparison of resins for fabricating provisional fixed restorations. Int J Prosthodont 1989;2:173-84.
- Balkenhol M, Ferger P, Mautner MC, Wöstmann B. Provisional crown and fixed partial denture materials: mechanical properties and degree of conversion. Dent Mater 2007;23:1574-83.
- 16. McLean JW. The failed restoration: causes of failure and how to prevent them. Int Dent J 1990;40:354-8.
- 17. Balkenhol M, Knapp M, Ferger P, Heun U, Wöstmann B. Correlation between polymerization shrinkage and marginal fit of temporary crowns. Dent Mater 2008;24:1575-84.
- Gougaloff R, Stalley FC. Immediate placement and provisionalization of a dental implant utilizing the CEREC 3 CAD/CAM Protocol: a clinical case report. J Calif Dent Assoc 2010;38:170-3, 176-7.
- Poticny DJ, Klim J. CAD/CAM in-office technology: innovations after 25 years for predictable, esthetic outcomes. J Am Dent Assoc 2010;141:5S-9S.
- Ozcan M, Nijhuis H, Valandro LF. Effect of various surface conditioning methods on the adhesion of dual-cure resin cement with MDP functional monomer to zirconia after thermal aging. Dent Mater J 2008;27:99-104.
- Pisani-Proenca J, Erhardt MC, Valandro LF, Gutierrez-Aceves G, Bolanos-Carmona MV, Del Castillo-Salmeron R, Bottino MA. Influence of ceramic surface conditioning and resin cements on microtensile bond strength to a glass ceramic. J Prosthet Dent 2006;96:412-7.
- 22. Alt V, Hannig M, Wöstmann B, Balkenhol M. Fracture strength of temporary fixed partial dentures: CAD/CAM versus directly fabricated restorations. Dent Mater 2011;27: 339-47.
- 23. Balkenhol M, Meyer M, Michel K, Ferger P, Wöstmann B. Effect of surface condition and storage time on the repairability of temporary crown and fixed partial denture materials. J Dent 2008;36:861-72.
- 24. Balkenhol M, Mautner MC, Ferger P, Wöstmann B. Mechanical properties of provisional crown and bridge materials: chemical-curing versus dual-curing systems. J Dent 2008;36: 15-20.
- 25. Kerby RE, Knobloch LA, Sharples S, Peregrina A. Mechanical properties of urethane and bis-acryl interim resin materials. J Prosthet Dent 2013;110:21-8.

- Knobloch LA, Kerby RE, Pulido T, Johnston WM. Relative fracture toughness of bis-acryl interim resin materials. J Prosthet Dent 2011;106:118-25.
- 27. Cheng CJ, Lin CL, Shan YF. Multifactorial analysis of variables influencing the fracture strength of repair joints for provisional restorative materials using the statistically based Taguchi method. J Dent Sci 2010;5:90-9.
- Nejatidanesh F, Momeni G, Savabi O. Flexural strength of interim resin materials for fixed prosthodontics. J Prosthodont 2009;18:507-11.
- 29. Lang R, Rosentritt M, Behr M, Handel G. Fracture resistance of PMMA and resin matrix composite-based interim FPD materials. Int J Prosthodont 2003;16:381-4.
- Pihut M, Wisniewska G, Majewski P, Gronkiewicz K, Majewski S. Measurement of occlusal forces in the therapy of functional disorders with the use of botulinum toxin type A. J Physiol Pharmacol 2009;60:113-6.
- Tortopidis D, Lyons MF, Baxendale RH, Gilmour WH. The variability of bite force measurement between sessions, in different positions within the dental arch. J Oral Rehabil 1998;25:681-6.