

Original article

eISSN 2234-0742 Genomics Inform 2020;18(4):e40 https://doi.org/10.5808/GI.2020.18.4.e40

Received: October 6, 2020 Revised: October 16, 2020 Accepted: October 16, 2020

*Corresponding author: E-mail: sunshin@catholic.ac.kr

**Corresponding author: E-mail: yejun@catholic.ac.kr

© 2020, Korea Genome Organization

C This is an open-access article distributed under the terms of the Creative Commons Attribution license (http://creativecommons. org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Development of reverse transcription loop-mediated isothermal amplification assays for point-of-care testing of avian influenza virus subtype H5 and H9

Songzi Zhang^{1,2,3,4}, Juyoun Shin⁵, Sun Shin^{2,3,4*}, Yeun–Jun Chung^{1,2,3,4**}

¹Department of Biomedicine & Health Sciences, Graduate School, The Catholic University of Korea, Seoul 06591, Korea

²Department of Microbiology, College of Medicine, The Catholic University of Korea, Seoul 06591, Korea

³Precision Medicine Research Center, College of Medicine, The Catholic University of Korea, Seoul 06591, Korea

⁴Integrated Research Center for Genome Polymorphism, College of Medicine, The Catholic University of Korea, Seoul 06591, Korea

⁵ConnectaGen, Hanam 12918, Korea

Avian influenza (AIV) outbreaks can induce fatal human pulmonary infections in addition to economic losses to the poultry industry. In this study, we aimed to develop a rapid and sensitive point-of-care AIV test using loop-mediated isothermal amplification (LAMP) technology. We designed three sets of reverse transcription LAMP (RT-LAMP) primers targeting the matrix (M) and hemagglutinin (HA) genes of the H5 and H9 subtypes. RT-LAMP targeting the universal M gene was designed to screen for the presence of AIV and RT-LAMP assays targeting H5-HA and H9-HA were designed to discriminate between the H5 and H9 subtypes. All three RT-LAMP assays showed specific amplification results without nonspecific reactions. In terms of sensitivity, the detection limits of our RT-LAMP assays were 100 to 1,000 RNA copies per reaction, which were 10 times more sensitive than the detection limits of the reference reverse-transcription polymerase chain reaction (RT-PCR) (1,000 to 10,000 RNA copies per reaction). The reaction time of our RT-LAMP assays was less than 30 min, which was approximately four times quicker than that of conventional RT-PCR. Altogether, these assays successfully detected the existence of AIV and discriminated between the H5 or H9 subtypes with higher sensitivity and less time than the conventional RT-PCR assay.

Keywords: avian influenza virus, H5 subtype, H9 subtype, RT-LAMP

Introduction

Influenza A viruses, which are classified into 18 hemagglutinin (HA) and 11 neuraminidase subtypes based on their antigenic properties, have been found to infect not only avian species, but also multiple mammalian hosts including humans. In particular, avian influenza viruses (AIVs) belong to the category of influenza A viruses, which induce human infections and cause economic losses annually [1-3].

Starting in late 2003, simultaneous outbreaks of highly pathogenic avian influenza virus

(HPAIV) in poultry occurred across diverse Asian countries, mainly caused by influenza H5N1 virus [2,3]. The outbreaks spread widely throughout the globe [4], and caused several human infections, including fatal cases [5]. In addition to HPAIVs, low-pathogenic avian influenza viruses (LPAIVs), which usually produce mild or no symptoms in birds, can induce pathogenic avian influenza by antigenic drift or shift, and can also induce human infection [2]. Indeed, the H9 subtype AIV, an LPAIV, was reported to be capable of human infection [6], and the World Health Organization warned that the H9N2 subtype could trigger a global influenza outbreak in humans, albeit with relatively lower pathogenicity than H5N1 or H7N9. Therefore, identification of HPAIVs and LPAIVs would be important to control avian influenza and potential human pandemic infections.

To identify AIVs at the early stage of disease, a rapid, specific, and sensitive detection method is required. Polymerase chain reaction (PCR) and reverse transcription PCR (RT-PCR) are the most commonly used tools, serving as a gold standard for the diagnosis of viral infections including influenza. However, these methods require specialized equipment and trained personnel for the entire experimental process [7], which makes these assays hard to use for point-of-care testing (POCT) in the field. In addition, ordinary PCR or RT-PCR requires about 2 hours to process and well-purified nucleic acids because the reaction is sensitive to PCR inhibitors.

Notomi et al. [8] developed a new nucleic acid amplification technology, loop-mediated isothermal amplification (LAMP), which can amplify target nucleic acids at a consistent temperature without changing the temperature for denaturation, annealing, and extension. LAMP only needs strand displacement by using the Bst DNA polymerase enzyme for nucleic acid amplification, so it can be amplified under isothermal conditions with a high amplification efficiency [8]. LAMP technology has a number of advantages for POCT [9-11]. For example, the time required for LAMP assays is much shorter than that for PCR, and LAMP assays can detect the target nucleic acid with high sensitivity without the equipment needed for thermal cycling. Furthermore, diverse detection methods can be applied, such as turbidity, colorimetric detection, and fluorescent dye incorporation.

Due to those advantages, LAMP technology has been widely applied to detect various pathogens and its applications have been advanced to include several different forms of assays, such as reverse transcription LAMP (RT-LAMP) and multiplex LAMP [12,13]. In particular, RT-LAMP assays have been found to be an efficient tool to detect RNA viruses such as the severe acute respiratory syndrome coronavirus 2, influenza, dengue, and West Nile viruses [14-22]. In RT-LAMP, reverse transcription of viral RNA and cDNA amplification can be processed in a single tube as one step under isothermal conditions. The straightforwardness of RT-LAMP is ideal for on-site field testing or diagnostics in a POCT setting where laboratory equipment is limited [13]. Together, RT-LAMP can be a simple and effective tool to detect and distinguish diverse AIV subtypes in the poultry field or at clinical quick-check desks [19].

In this study, we developed RT-LAMP assays that can detect a universal target for AIV and can discriminate between the H5 and H9 subtypes of AIV. We validated the sensitivity, specificity, and processing time of our assay.

Methods

Primer design for RT-LAMP assays

The nucleotide sequences of the matrix (M) gene and the HA genes from the H5 and H9 subtypes (H5-HA and H9-HA, respectively) were retrieved from the Influenza Research Database (IRD) and Global Initiative on Sharing All Influenza Data (GISAID), and aligned using DNAStar (Lasergene, Madison, WI, USA). The RT-LAMP assay primers were designed based on M, H5-HA and H9-HA sequence alignments by using the software PrimerExplorer version 5 (http://primerexplorer.jp/lampvSe/). The feasibility of all sets of primers was then subsequently validated using the BLAST program.

Preparation of target RNAs

Full-length fragments of the AIV M, H5-HA, and H9-HA genes with the T7 promoter sequence (TAATACGACTCACTATAGG-GAGA) were chemically synthesized and cloned into a pTwist Amp High Copy plasmid (Twist Bioscience, South San Francisco, CA, USA). The RNA was transcribed by the universal M13F and M13R-pUC primers using the T7 RiboMAX Express Large Scale RNA Production System (Promega, Madison, WI, USA). The transcribed RNA was 10-fold serially diluted from 10° copies/µL to 10 copies/µL, and used as a template for RT-LAMP and RT-PCR.

RT-LAMP

The RT-LAMP reaction was carried out as described elsewhere [18-20]. In total, 25 μ L of a mixture containing 0.2 μ M of each outer primer (F3 and B3), 1.6 μ M of each inner primer (FIP and BIP), 0.8 μ M of each loop primer (LF and LB), 8 U of the Bst 2.0 DNA polymerase (New England Biolabs, Hitchin, UK), 2 U of AMV reverse transcriptase (New England Biolabs), 8 mM of MgSO₄ (New England Biolabs), 1.4 mM of each dNTPs (Thermo Fisher Scientific, Waltham, MA, USA), 1 × isothermal amplification buffer (New

England Biolabs) and 0.4 M N-methylformamide (NMF) and isobutylamide (IBA) was prepared. For the real-time assay, 1 U of SYTO9 stain (Thermo Fisher Scientific) was added. One microliter of transcribed RNA was added to the respective tube. The reaction was carried out at 60°C for the H5-HA gene or at 68°C for the M gene and the H9-HA gene in a thermal cycler or CFX96 Touch Real-Time PCR Detection System (Bio-Rad Laboratories, Hercules, CA, USA). The reaction was terminated by heating at 80°C for 10 min. The fluorescence curve was captured in real time, and the result was illustrated as a graph on the monitor of the real-time system, verifying the amplification. RT-LAMP products were then evaluated by electrophoresis using 1.5% agarose gels to ensure the products of amplification reaction.

RT-PCR

One-step RT-PCR was performed with the AIV Multi-tube RT-PCR Kit (iNtRON Biotechnology, Seongnam, Korea), which was approved by the Animal and Plant Quarantine Agency of Republic of Korea, according to the manufacturer's instructions. The reaction cycling conditions were as follows: 30 min of RT at 45°C, RT inactivation and polymerase activation for 5 min at 94°C, 40 cycles at 94°C for 30 s, 55°C for 60 s and 72°C for 60 s, and final extension for 5 min at 72°C. The RT-PCR products were then evaluated by electrophoresis using 1.5% agarose gel. The expected product sizes were 378 bp for the M gene, 311 bp for the H5-HA gene, and 252 bp for the H9-HA gene, respectively.

Results

Primer design and optimization of RT-LAMP assays

For RT-LAMP, primers were designed for the M, H5-HA, and H9-HA genes. The six primer sets for each target gene comprised two outer primers (F3 and B3), two inner primers (FIP and BIP), and two loop primers (LF and LB). In this study, in vitro-transcribed full-length RNA fragments of the M, H5-HA, and H9-HA genes were used as templates. Two to three different primer sets for each gene were tested and the sets that demonstrated the best amplification performance were selected for further analysis (Supplementary Fig. 1). The sequence information of the final RT-LAMP primer sets is available in Table 1. All selected primers in this study were designed within a conserved region of AIVs isolated from diverse poultry (Fig. 1). To minimize nonspecific amplification, we tested four different combinations of NMF and IBA in the reaction mixture (0.2 M NMF + IBA, 0.4 M NMF + IBA, 0.6 M NMF + IBA, 0.8 M NMF + IBA) using the H9-HA RT-LAMP set. When we applied a specific target template (H9-HA), all four reaction conditions showed significant amplification signals within 20 min regardless of the combinations of the additives (Fig. 2). When we applied a negative control template, the reaction mixtures with 0.4, 0.6, or 0.8 M NMF + IBA did not show any nonspe-

 Table 1. Primer sets designed for the detection of the M gene and H5- and H9-HA genes of AIVs

Gene	Primer	Length (bp)	Sequence (5'-3')
Μ	F3	19	GCATCGGTCTCACAGACAG
	B3	19	ACTGGAGCTAGGGTGAGTT
	FIP ^a	46	CAGCCATCTGCTCCATAGCCTTTTTTCCACCAACCCACTAATCAGG
	BIP ^a	42	CAGCGGAAGCCATGGAGGTTTTTCCAATGGTCCTCATCGCC
	LPF	19	GGCCAGCACCATTCTGTTT
	LPB	18	AGGCTAGGCAGATGGTGC
H5-HA	F3	23	GCATACAAAATTGTCAAGAAAGG
	B3	19	ACTATTTCTGAGTCCAGTC
	FIP ^a	52	CGCMCCTATTGGAGTTTGACATTTTTCTCAACAATTATGAAAAGTGA
	BIP ^a	52	TCTAGYATGCCRTTCCACAATATACTTTTGCAAGGACTAATTTGTTTGATTT
	LF	20	GTGTTGCAGTGGCCATACTC
	LB	21	ACCCTCTCACCATCGGGGAAT
H9-HA	F3	19	GAGAATCCTGAAGACCGAT
	B3	19	CCCTCTACCTGATGTAGCA
	FIP ^a	48	GTGGAATGGCAATGTCGTATTCAAATTTTTTAAATAGTGGCAACTGCG
	BIP ^a	46	AGCAAGTATGCATTTGGGAACTTTTTGCACATTTCTTAGACCAACT
	LF	20	CCTTTCTCAGTTTGGCACTG
	LB	25	GTTGGAGTGAAGAGTCTCAAACTGG

M, matrix; HA, hemagglutinin; AIV, avian influenza virus; RT-LAMP, reverse transcription loop-mediated isothermal amplification. ^aEach inner primer (FIP and BIP) of RT-LAMP had two binding regions (F1c + F2 and B1c + B2, respectively) connected by a TTTT spacer.

A	F3	F2	LF	
	GCATCGGTCTCACAGACAG	CCACCAACCCACTAATCAGG	AAACAGAATGGTGCTGGCC	AAGGCTA
A/duck/Wuhan/WHYF05/2014	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	********************	*******
A/duck/Bangladesn/24075/2015 A/teal/Egypt/MB-D-698OP/2016	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	***********************************	********
A/wild waterfowl/Korea/F7-18/2018	* * * * * * * * * * * * * * * * * * * *	*T * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	*******
A/domestic duck/Georgia/1/2016	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	***********************	******
A/European teal/Chany/63/2011 A/duck/Viet Nam/HU10-48/2018	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	***********************************	*********
A/mallard duct/Netherlands/60/2015	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	******
A/duck/Mongolia/914/2018	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	*******************	*******
A/duck/Mongolia/786/2017	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	**********************	*******
F1	B1	LB	B2	
1GGAGCAGA1GGC1GC.	AGCGGAAGCCA I GGAGG I I/	AGGCTAGGCAGATGGTGC-GGC	GATGAGGACCATTGGAACTCACC	CTAGCTCCAGT
* * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * *
* * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * *
* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * *
* * * * * * * * * * * * * * * * * *	* * *A * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * *
* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * A * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * *
* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * *
* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * * *	****G******G******	* * * * * * * * * * *
В	F3	F2	LF	F1
A /dwat/Eastana China (00404/0044	GCATACAAAATTGTCAAGAAAGG	CTCAACAATTATGAAAAGTGA	GAGTATGGCCACTGCAACAC	AATGTCAAACT
A/chicken/Montana/15-010559-1/2015	*********************	* * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	* * * * * * * * * * * *
A/breeder chicken/Korea/H250/2014	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	* * * * * * * * * * *
A/chicken/Zhejiang/209/2016	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	* * * * * * * * * * *
A/breeder chicken/Korea/H122/2014	*********************	* * * * * * * * * * * * * * * * * *	* * # * * C * * * * * * * * * * * * * *	* * * * * * * * * * * *
A/chicken/Netherlands/emc-3/2014	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * *	* * * * * * * * * * C
A/muscovy duck/Viet Nam/HN-1697/2014	* * * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * *	* * A * * * * * * * * * * * * * * *	* * * * * * * * * * *
A/goose/Taiwan/01006/2015	* * * * * * * * * * * * * * * * * * * *	****	* * * * * * * * * * * * * * * * *	* * * * * * * * * * *
A/WHOODEL SWAII/ LOCHIGI/ 1/2017	**********A********	* * * * * * * * * * * * * * * * * *	* *A * * * * * * * * * * * * * * *	* * * * * * * * * *
Avwhooper swan/ rochig// 1/2017	************A***********	****	* *A * * * * * * * * * * * * * * *	· · · · · · · · · · · · · · ·
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	LB CACAATATACACCCTCTCACCAT	**********************************	* A * * * * * * * * * * * * * * * *	
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	LB CCACAATATACACCCTCTCACCAT	AAATCAAA	B2 L CAAATTAGTCCTTGCGAACTGGAC G	- B TCAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	LB CACAATATACACCCTCTCACCAT		B2 L CAAATTAGTCCTTGC/GACTGGAC G	- B TCAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	A LB CACAATATACACCCTCTCACCAT	CGGGGAATAAATCAAA - G. - G.	B2 L CAAATTAGTCCTTGC/BACTGGAC G G G	-B TCAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	A LB CACAATATACACCCTCTCACCAT	2GGGGAATAAATCAAA - G - G 	B2 L CAAATTAGTCCTTGCBACTGGAC G G G	-B TCAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	A LB CACAATATACACCCTCTCACCATC	2GGGGAATAAATCAAA 	B2 L CAAATTAGTCCTTGCBACTGGAC G G G G	-B TCAGAAATAGT
ECCAATAGGKGCGTCTAGYATGCCRTTC	А	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G	-B TCAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	A LB CACAATATACACCCTCTCACCATO	CGGGGAATAAATCAAA 	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	-B TCAGAAATAGT
AWNOOPEN SWAIN TOCINGUIN2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC	LB CACAATATACACCCTCTCACCAT	CGGGGAATAAATCAAA - G - G - G - G - G - G - G	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G C G G C G G C G C C G	-B TCAGAAATAGT
AWNOUPER SWAIN TOCINGUIN2017	LB CACAATATACACCCTCTCACCATO	CGGGGAATAAATCCAAA 	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	-B ICAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	EB CACAATATACACCCTCTCACCATO	CGGGGAAT AAATCAAA 	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G C G C G C	.B ICAGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC	В САСААТАТАСАСССТСТСАССАТ 	CGGGGAAT	B2 L CAAATTAGTCCTTGC/SACTGGAC G G G G G G G G G C G C G C G C C C C	-B ICAGAAATAGT F1 TACGACATTGC
Aveniouper swann rocing: h2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC TT CCAATAGKGCGTCTAGYATGCCRTTC Avenioused for the state of t	EB ICACAATATACACCCTTCTCACCATO IIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIII	CGGGGAAT	B2 LF CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G	-B ICAGAAATAGT F1 TACGACATTGC
Advinooper swanrocing/ h2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC 	В сасаататасассстстсассато т. 	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	-B ICAGAAATAGT F1 TACGACATTGC
Advinoupler swall rocing/ h2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC 	EB ECACAATATACACCCTCTCACCAT T. 	CGGGGAATAAATCAAA 	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	В ГСАБАААТАБТ F1 TACGACATTGC
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC A/chicken/Korea/C3/2009 A/chicken/Korea/SH0915/2009 A/chicken/KOP15/2009	В Сасаататасассстстсассато 	CGGGGAATAAATCCAAA 	B2 L CAAATTAGTCCTTGC/BACTGGAC G G G G G G G G G G G G G G G G G G	-B TCAGAAATAGT F1 TACGACATTGC
Avinooper swannooling/ h2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC 	В ССАСААТАТАСАСССТСТСАССАТ 	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	-B ICAGAAATAGT F1 TACGACATTGC
Avmiouper swannochigi h2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCATAGGKGCGTCTAGYATGCCRTTC Avchicken/Korea/C3/2009 Avchicken/Korea/C3/2009 Avchicken/Korea/KNUWSJ09/2009 Avchicken/Korea/KNU915/2009 Avchicken/Korea/KU1627/2007 Avchicken/Korea/KU1627/2007 Avchicken/Korea/KU10-194/2007 Avchicken/Korea/KU10-194/2007	В сасаататасассстстсассата т.	CGGGGAAT	B2 LF GCCAAACTGAGAGAGAGAGAGAGAGAGAGAGAGAGAGAGA	-B ICAGAAATAGT F1 TACGACATTGC
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC 	EB ICACAATATACACCCTTCTCACCATO F3 GAGAATCCTGAAGACCGAIIITAA	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	F1 TACGACATTGC
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC 	EB ICACAATATACACCCTCTCACCATO F3 GAGAATCCTGAAGACCGATIITAA	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGACC G G G G G G G G G G G G G G G G G	F1 TACGACATTGC
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC 	В САСААТАТАСАСССТСТСАССАТО 	CGGGGAATAAATCAAA 	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G	F1 TACGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC 	В ССАСААТАТАСАСССТСТСАССАТ	CGGGGAATAAATCAAA 	В2 Ц САААТТАGTCCTTGC/GACTGGACC G G G G G G G G G G G G G G G G G	E TCAGAAATAGT F1 TACGACATTGC T
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC 	В САСААТАТАСАСССТСТСАССАТА 	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G	F1 TACGACATTGC T.CCGACATTGC
Avmiouper swannocing/ h2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC Archicken/Korea/C3/2009 Archicken/Korea/KNUWSJ09/2009 Archicken/Korea/KNUSJ09/2009 Archicken/Korea/SH0915/2009 Archicken/Korea/SH0915/2009 Archicken/Korea/SH0915/2009 Archicken/Korea/SH0915/2009 Archicken/Korea/SH0915/2007 Archicken/Korea/SH0915/2007 Archicken/Korea/KU10-194/2007 Archicken/Korea/At5/2006 Archicken/Korea/At5/2006 Archicken/Korea/At5/2006 Archicken/Korea/At5/2007 B1 CATTCCACAGCAAGTATGCAT	EB ICACAATATACACCCTCTCACCATO TT F3 GAGAATCCTGAAGACCGATITTAA GAGAATCCTGAAGACCGATITTAA	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G G G G G G	F1 TACGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC A/chicken/Korea/C3/2009 A/chicken/Korea/KNUWSJ09/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/SH0915/2007 A/chicken/Korea/A08/2007 B1 CATTCCACAGCAAGTATGCAT	F3 GAGAATCCTGAGAGACCGATIITAA	CGGGGAAT	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G	F1 TACGAAATAGT F1 TACGACATTGC
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC A/chicken/Korea/C3/2009 A/chicken/Korea/KNUWSJ09/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2009 A/chicken/Korea/SH0915/2007 A/chic	EB ICACAATATACACCCTCTCACCATO F3 GAGAATCCTGAAGACCGATIITAA GAGAATCCTGAAGACCGATIITAA	CGGGGAATAAATCAAA 	В2 L САААТТАGTCCTTGC/GACTGGACC G	F1 TACGAAATAGT
B1 CCAATAGGKGCGTCTAGYATGCCRTTC CCAATAGGKGCGTCTAGYATGCCRTTC 	F3 GAGAATCCTGAAGACCGATIITAA	CGGGGAATAAATCAAA 	B2 LF GCCAAACTGAGAAAGGTTTGAA B2 LF GCCAAACTGAGAAAGGTTTGAA B2 L B2 L B2 L B2 L B2 L B2 L B2 L B3 L	B TCAGAAATAGT F1 TACGACATTGC T B AGGTAGAGGG
Advindoper swann rochigi //2017 B1 CCAATAGGKGCGTCTAGYATGCCRTTC Advincken/Korea/C3/2009 Advincken/Korea/KNUWSJ09/2009 Advincken/Korea/KNU95/2009 Advincken/Korea/KU909/2007 Advincken/Korea/KU10-194/2007 Advi	EB LB LB LB LB LB LB TIGGGAATTATACACCCTOTCACCATO T T T T T T T T T T T T T	CGGGGAATAAATCAAA 	B2 L CAAATTAGTCCTTGC/GACTGGAC G G G G G G G G G G G G G	EB ICAGAAATAGT F1 TACGACATTGC

Fig. 1. Primer design for the reverse transcription loop-mediated isothermal amplification assays. Multiple sequence alignments of the M (A), H5-specific (B), and H9-specific (C) HA genes with 10 nucleotide sequences from BLAST. The F3, F2, LF, F1, B1, LB, B2, and B3 regions are indicated above the sequences. M, matrix; HA, hemagglutinin.

cific signals during the whole process (60 min); however, the reaction containing 0.2 M NMF + IBA showed a subtle nonspecific signal after 50 min (Fig. 2). Therefore, we decided to add 0.4 M NMF + IBA for all three RT-LAMP reactions.

Specificity of the RT-LAMP assays

To evaluate the specificity of the three RT-LAMP assays, we ap-

plied matched and non-matched template RNA samples including a negative control (no template). The RT-LAMP assay for the universal M target showed an amplification product with the M RNA template, but no amplification product was detected with for the HA RNA templates (H5 and H9) or the negative control template (Fig. 3A). The RT-LAMP assays for H5-HA and H9-HA showed target-specific amplification products without any



Fig. 2. Optimization of the combinations of additives (NMF and IBA) in the reverse transcription loop-mediated isothermal amplification reaction mixture. (A) 0.2 M NMF + IBA. (B) 0.4 M NMF + IBA. (C) 0.6 M NMF + IBA. (D) 0.8 M NMF + IBA. The x-axis represents the time for RT-LAMP reaction; the y-axis represents the relative fluorescence signal. NMF, N-methylformamide; IBA, isobutylamide; NTC, negative control; RT-LAMP, reverse transcription loop-mediated isothermal amplification.



Fig. 3. Specificity of the RT-LAMP assays. (A) The RT-LAMP products of the RT-LAMP assay for M gene were electrophoresed with 1.5% agarose gel: lane 1, negative control; lane 2, synthesized M RNA template; lane 3, synthesized H5-HA RNA template; lane 4 synthesized H9-HA RNA template; lane M, 100-bp DNA marker. (B, C) RT-LAMP products for H5-specific and H9-specific HA genes were electrophoresed with 1.5% agarose gel, respectively: lane 1, negative control; lane 2, synthesized H5-HA RNA template; lane 3, synthesized H9-HA RNA template; lane M, 100-bp DNA marker. RT-LAMP, reverse transcription loop-mediated isothermal amplification; M, matrix; HA, hemagglutinin.

cross-reaction with non-matched templates or the negative control (Fig. 3B and 3C).

Sensitivity of the RT-LAMP assays

To evaluate the sensitivity of the assays, we examined the detection limits of the three RT-LAMP assays. To do so, the template RNAs were 10-fold serially diluted (ranging from 10 to 10⁵ RNA copies) and applied for each real-time RT-LAMP assay. The detection limit of the M gene and the H5-HA RT-LAMP assays was found to be 100 copies/reaction, while the detection limit of the H9-HA RT-LAMP assay was 1,000 copies/reaction (Fig. 4A-4C). For objective validation of the sensitivity of our RT-LAMP assay, we performed conventional RT-PCR with a certified commercial AIV detection kit (AIV Multi-tube RT-PCR Kit, iNtRON Biotechnology) using the same template RNAs used for the RT-LAMP assays and compared the results. Our RT-LAMP assay was 10 times more sensitive than the reference RT-PCR kit. The detection limit of the RT-PCR assays for the M and H5-HA genes was 1,000 copies/reaction, and that for the H9-HA gene was 10,000 copies/reaction (Fig. 4D). Regarding the reaction time, the fluorescence signal of target-specific amplification appeared within 20 min in all three RT-LAMP assays when 10⁵ RNA copies were used. In the RT-

LAMP assay for the M gene, even with 100 copies of the RNA template (the limit of detection), the amplification signal appeared within 20 min.

Discussion

Recently emerging viral infectious diseases, including AIV, are increasing, posing a major threat to both public health and poultry farming. POCT is a new concept of laboratory testing that enables testing to be performed where an infection occurs without transporting the samples to central clinical laboratories [7]. LAMP technology is ideal for POCT due to its ease of performance without the need for sophisticated equipment or experts to operate it [13]. In particular, RT-LAMP can be useful to identify RNA viruses at remote locations where laboratory equipment is limited [12,13]. In this study, we aimed to developed rapid and highly sensitive assays for identifying AIVs through POCT in the context of monitoring AIV epidemic outbreaks. Using LAMP technology, we developed RT-LAMP assays that can detect the presence of AIV and distinguish the H5 and H9 subtypes within 20-30 minutes. Furthermore, all three RT-LAMP assays (M, H5, and H9) were approximately 10-fold more sensitive than the approved ref-



Fig. 4. Sensitivity of RT-LAMP and RT-PCR. The template RNAs, M (A), H5-HA (B), and H9-HA (C), were 10-fold serially diluted (ranging from 10 to 10⁵ RNA copies) and applied for each real-time RT-LAMP assay. The x-axis represents the time for the RT-LAMP reaction; the y-axis represents the relative fluorescent signal. (D) The same template RNAs were applied for conventional RT-PCR using a certified commercial AIV detection kit (AIV Multi-tube RT-PCR Kit, iNtRON Biotechnology) and the products were electrophoresed with 1.5% agarose gel. RT-LAMP, reverse transcription loop-mediated isothermal amplification; RT-PCR, reverse transcription polymerase chain reaction; M, matrix; HA, hemagglutinin; AIV, avian influenza virus; NTC, negative control.

erence RT-PCR assay.

We first designed an AIV universal RT-LAMP assay by targeting the M gene, which has a common sequence across AIV subtypes. Therefore, a positive M RT-LAMP assay result indicated the presence of AIV in the sample. We also designed RT-LAMP primer sets specifically targeting H5-HA and H9-HA, as these are typical subtypes of HPAIVs and LPAIVs, respectively.

All RT-LAMP primers were designed within a conserved region of AIVs isolated from diverse poultry across the world, suggesting that our RT-LAMP assays can be suitable for determining the presence of HS-HA and H9-HA regardless of the poultry or location. After designing the RT-LAMP primers, we used *in vitro*-transcribed full-length RNA fragments of the M, H5-HA, and H9-HA genes for experimental validation, because field-isolated AIV samples were not available due to their high pathogenicity.

Due to the very high level of amplification efficiency of the LAMP, false detection is a major concern with this technology. In this study, we used additives (NMF and IBA) to minimize non-specific amplification [23] and checked the real-time amplification signal by adding a fluorescent reagent. Finally, we selected 0.4 M NMF + IBA because this condition showed an efficient target-specific amplification signal without any nonspecific amplification signal. Moreover, none of the three RT-LAMP assays showed any nonspecific reactions with non-matched template or negative controls. These results suggest that our RT-LAMP assays are specific enough to identify AIV subtypes.

When we checked the sensitivity, the detection limits of our RT-LAMP assays were 100 to 1,000 RNA copies per reaction, which were 10 times more sensitive than the approved reference RT-PCR assays (1000 to 10,000 RNA copies per reaction). These results are consistent with previous studies reporting that LAMP assays had higher sensitivity than PCR-based assays [18,19]. In particular, the reaction time of our RT-LAMP assay was less than 30 minutes, which was approximately four times quicker than that of the conventional RT-PCR assay. These results suggest that our RT-LAMP assays can identify AIV subtypes 10 times more sensitively and four times more quickly than the conventional RT-PCR assays.

In summary, the RT-LAMP assays targeting the M, H5-HA, and H9-HA genes developed in this study demonstrated a high level of sensitivity and specificity. The assays could successfully detect the existence of AIV and discriminate between the H5 and H9 subtypes with higher sensitivity and less time than the conventional RT-PCR assays. This method could be a useful POCT tool for the rapid identification of AIV infections in the field.

ORCID

Songzi Zhang: https://orcid.org/0000-0003-4641-499X Juyoun Shin: https://orcid.org/0000-0001-9207-7053 Sun Shin: https://orcid.org/0000-0002-2146-4336 Yeun-Jun Chung: https://orcid.org/0000-0002-6943-5948

Authors' Contribution

Conceptualization: SS, JS, YJC. Data curation: SZ, SS, YJC. Formal analysis: SZ, SS. Funding acquisition: YJC. Methodology: SZ, JS, SS. Writing - original draft: SZ, SS. Writing - review & editing: YJC.

Conflicts of Interest

No potential conflict of interest relevant to this article was reported.

Acknowledgments

This work was supported by a grant from the National Research Foundation of Korea (2017R1E1A1A01074913 and 2017M3C9 A6047615). We thank KREONET (Korea Research Environment Open NETwork) and KISTI (Korea Institute of Science and Technology Information) for allowing us to use their network infrastructure.

Supplementary Materials

Supplementary data can be found with this article online at http://www.genominfo.org.

References

- 1. Webster RG, Bean WJ, Gorman OT, Chambers TM, Kawaoka Y. Evolution and ecology of influenza A viruses. Microbiol Rev 1992;56:152-179.
- 2. Kim SM, Kim YI, Pascua PN, Choi YK. Avian influenza A viruses: evolution and zoonotic infection. Semin Respir Crit Care Med 2016;37:501-511.
- Medina RA, Garcia-Sastre A. Influenza A viruses: new research developments. Nat Rev Microbiol 2011;9:590-603.
- 4. Office International des Epizooties (OIE). Update on avian influenza in animals (types H5 and H7). Paris: World Organisation for Animal Health, 2020. Accessed 2020 Jul 16. Available from: https://www.oie.int/en/animal-health-in-the-world/update-on-avian-influenza.

- 5. World Health Organization. Cumulative number of confirmed human cases of avian influenza A (H5N1) reported to WHO, 2003-2020. Geneva: World Health Organization, 2020. Accessed on 2020 May 8. Available from: https://www.who.int/influenza/ human_animal_interface/2020_MAY_tableH5N1.pdf?ua = 1. pdf.
- 6. Yu X, Jin T, Cui Y, Pu X, Li J, Xu J, et al. Influenza H7N9 and H9N2 viruses: coexistence in poultry linked to human H7N9 infection and genome characteristics. J Virol 2014;88:3423-3431.
- Niemz A, Ferguson TM, Boyle DS. Point-of-care nucleic acid testing for infectious diseases. Trends Biotechnol 2011;29:240-250.
- Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K, Amino N, et al. Loop-mediated isothermal amplification of DNA. Nucleic Acids Res 2000;28:E63.
- **9.** Kashir J, Yaqinuddin A. Loop mediated isothermal amplification (LAMP) assays as a rapid diagnostic for COVID-19. Med Hypotheses 2020;141:109786.
- Shirato K, Yano T, Senba S, Akachi S, Kobayashi T, Nishinaka T, et al. Detection of Middle East respiratory syndrome coronavirus using reverse transcription loop-mediated isothermal amplification (RT-LAMP). Virol J 2014;11:139.
- Miyamoto S, Sano S, Takahashi K, Jikihara T. Method for colorimetric detection of double-stranded nucleic acid using leuco triphenylmethane dyes. Anal Biochem 2015;473:28-33.
- 12. Mori Y, Notomi T. Loop-mediated isothermal amplification (LAMP): a rapid, accurate, and cost-effective diagnostic method for infectious diseases. J Infect Chemother 2009;15:62-69.
- Wong YP, Othman S, Lau YL, Radu S, Chee HY. Loop-mediated isothermal amplification (LAMP): a versatile technique for detection of micro-organisms. J Appl Microbiol 2018;124:626-643.
- 14. Kitagawa Y, Orihara Y, Kawamura R, Imai K, Sakai J, Tarumoto N, et al. Evaluation of rapid diagnosis of novel coronavirus disease (COVID-19) using loop-mediated isothermal amplification. J

Clin Virol 2020;129:104446.

- 15. Yan C, Cui J, Huang L, Du B, Chen L, Xue G, et al. Rapid and visual detection of 2019 novel coronavirus (SARS-CoV-2) by a reverse transcription loop-mediated isothermal amplification assay. Clin Microbiol Infect 2020;26:773-779.
- 16. Park GS, Ku K, Baek SH, Kim SJ, Kim SI, Kim BT, et al. Development of reverse transcription loop-mediated isothermal amplification assays targeting severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). J Mol Diagn 2020;22:729-735.
- Huang WE, Lim B, Hsu CC, Xiong D, Wu W, Yu Y, et al. RT-LAMP for rapid diagnosis of coronavirus SARS-CoV-2. Microb Biotechnol 2020;13:950-961.
- Poon LL, Leung CS, Chan KH, Lee JH, Yuen KY, Guan Y, et al. Detection of human influenza A viruses by loop-mediated isothermal amplification. J Clin Microbiol 2005;43:427-430.
- 19. Imai M, Ninomiya A, Minekawa H, Notomi T, Ishizaki T, Tashiro M, et al. Development of H5-RT-LAMP (loop-mediated isothermal amplification) system for rapid diagnosis of H5 avian influenza virus infection. Vaccine 2006;24:6679-6682.
- 20. Jang WS, Lim DH, Nam J, Mihn DC, Sung HW, Lim CS, et al. Development of a multiplex isothermal amplification molecular diagnosis method for on-site diagnosis of influenza. PLoS One 2020;15:e0238615.
- 21. Teoh BT, Sam SS, Tan KK, Johari J, Danlami MB, Hooi PS, et al. Detection of dengue viruses using reverse transcription-loop-mediated isothermal amplification. BMC Infect Dis 2013;13:387.
- 22. Parida M, Posadas G, Inoue S, Hasebe F, Morita K. Real-time reverse transcription loop-mediated isothermal amplification for rapid detection of West Nile virus. J Clin Microbiol 2004;42:257-263.
- 23. Tanner N, Evans TC. Compositions and methods for reducing background DNA amplification. US patent 9,546,358B2. 2017 Jan 17.