Photostability testing of IY-81149, a new Proton Pump Inhibitor

Kim SGO, Cheong JY, Jung WT, Kim DY and Cho KD

IlYang Central Research Institute, IlYang Pharmaceutical Co. LTD.

The study of the degradation of new drug substances under light is an integral part of stress testing because the photodegradation products might appear in the formulation during storage. These actinometric studies are used to identify precautionary measures and specify storage conditions needed in manufacturing or in formulation of the product. Also, these investigations provide valuable information for submission in applications for marketing authorizations for new active substances. Photostability testing of IY-81149, a new PPI substance, exposed to two artificial daylight lamps which have both UV and visible outputs (Dulux EL) and only visible light (Halogen Haloline), under the influence of container types such as amber, clear glass and quartz was studied. Both lamps provided overall illumination of not less than 1.2 million lux hours. We examined for any changes in color of IY-81149 substance and for assay and degradants arose from photochemical degradation processes. The result shows that the order of photoprotection of its container types were amber, quartz and clear glass. IY-81149 substance in amber was found to be photostable because it had no changes in color and was within justified limits ranged not more than 0.1% for each degradants and not more than 1% for total amount of degradants by the ICH stability and impurity guidelines. But in case of clear glass, this substance was photounstable manifesting that its contents decreased less than 99% after 3 weeks with above exposure to both light sources. We concluded the best storage container for IY-81149 was the amber in order to protect from photodegradation if the storage period was more than 3 weeks.

[PD4-5] [04/21/2000 (Fri) 14:50 - 15:50 / [1st Fl, Bldg 3]]

A Study on Pretreatment Method of Grains for Pesticide Residues Analysis (1)

Hwang YSo, Kim YS, Lee JM, Hwang KH, Kim JG, Kim MH, Seo SH

Kyung-Dong Agrochemical Analysis Team °, Seoul Metropolitan Government Institute of Health and Environment. College of Pharmacy, Kyung-Hee University*

This study was carried out to investigate optimal pretreatment condition for the pesticide residues in grains, including rice, barley, sticky rice etc. We bought 4 kinds of pesticide free grains and made high concentration contaminated grains (Metalaxyl, Parathion) artificially in lab. Diazinon was selected internal standard for GC-NPD and we confirmed pesticide residues spectrum by GC-mass. Accumulated pesticide toxicity was great dangerous, so need a rapid, accurate analysis method for grains. To get a better result, we used variation of time and wetting temperature, and recovery was not bad. According to kinds of grains, time, temperature and pesticide, different results were shown. So we have to more investigate a lot of kinds grains and pesticide for civil health and generation was used to consume great amounts grain daily

[PD4-6] [04/21/2000 (Fri) 14:50 - 15:50 / [1st Fl, Bldg 3]]

Chiral separation of β 2-agonists by Capillary electrophoresis using β -cyclodextrins as chiral selectors in a polyethylene glycol gel.

Kim KH1, Kim JH1, Kim HJ1, Jeun EYO1, Mar WC2, Kang JS3, Hong SP4

1College of Pharmacy, Kangwon National University, 2Natural Product Research Institute, Seoul National University, 3College of Pharmacy, Chungnam National University, 4College of Pharmacy,