Oral Presentation IV-1



Biologic responses of anodized titanium implants under different current voltage

최 정 원*, 허 성 주, 곽 재 영, 장 익 태, 김 현 이 ┃ 서울대학교 치과대학 보철학교실/공과대학 소재공학부

The oxide layer that covers titanium is very stable and seems to result in an excellent biocompatibility and successful osseointegration. The purpose of this study is to investigate the effect of high anodic oxidation voltage on the surface characteristics of titanium implants and biologic response on rabbit tibiae. Bone tissue responses were evaluated by removal torque test and histomorphometric analysis.

Screw shaped implants with microthreads (Warantec Co.) were made of commercially pure titanium(GradeII). We prepared anodized implant under 300V(groupI), 400V(groupII), 500V(groupIII) 550V(groupIV).

The surface characteristics of specimens were inspected in 3 categories: surface morphology, surface roughness and the thickness of oxide layer.

Screw shaped implants were installed in rabbit tibiae and the removal torque values were measured and histomorphometric analysis was done after 1 and 3months healing period.

From this study, we concluded that

- 1. As anodic oxidation voltage increased above 300V, oxide layer thickness increased rapidly and pore size increased.
- 2. The Ra values of group III and IV samples were higher than that of groupl and II samples. (P(0.05)
- 3. In removal torque test, groupIII after 1month healing period showed higher values than groupl and II in statistically significant level. (P(0.05)
- 4. In histomorphometric analysis groupIII and IV showed more BIC of total implant surface than groupl and II after 3months healing period. (P(0.05)

Acknowledgement: This work was supported by a grant from the Korea Health 21 R&D Project, Ministry of Health & Welfare, Republic of Korea (02-PJ3-PG6-EV11-0002)

Key word: anodic oxidation, voltage, histomorphometry, removal torque, titanium implant