

Application note

elSSN 2234-0742 Genomics Inform 2023;21(1):e13 https://doi.org/10.5808/gi.22076

Received: December 15, 2022 Revised: January 6, 2023 Accepted: January 6, 2023

*Corresponding author: E-mail: kyudong.han@gmail.com

Jinuk Jeong and Yunseok Oh contributed equally to this work.

© 2023 Korea Genome Organization

© This is an open-access article distributed under the terms of the Creative Commons Attribution license (http:// creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Effective microbial molecular diagnosis of periodontitis-related pathogen *Porphyromonas gingivalis* from salivary samples using *rgpA* gene

Jinuk Jeong¹, Yunseok Oh¹, Junhyeon Jeon², Dong-Heon Baek³, Dong Hee Kim⁴, Kornsorn Srikulnath⁵, Kyudong Han^{1,6,7,8}*

¹Department of Bioconvergence Engineering, Dankook University, Yongin 16890, Korea ²Department of Biological Sciences, Dankook University, Cheonan 31116, Korea

³Department of Oral Microbiology and Immunology, College of Dentistry, Dankook University, Cheonan 31116, Korea

⁴Department of Anesthesiology and Pain Management, Dankook University Hospital, Cheonan 31116, Korea

⁵Animal Genomics and Bioresource Research Unit (AGB Research Unit), Faculty of Science, Kasetsart University, Bangkok 10900, Thailand

⁶Center for Bio Medical Engineering Core Facility, Dankook University, Cheonan 31116, Korea ⁷Department of Microbiology, College of Science & Technology, Dankook University, Cheonan 31116, Korea

⁸HuNbiome Co., Ltd, R&D Center, Seoul 08503, Korea

Importance of accurate molecular diagnosis and guantification of particular disease-related pathogenic microorganisms is highlighted as an introductory step to prevent and care for diseases. In this study, we designed a primer/probe set for quantitative real-time polymerase chain reaction (gRT-PCR) targeting rgpA gene, known as the specific virulence factor of periodontitis-related pathogenic bacteria 'Porphyromonas gingivalis', and evaluated its diagnostic efficiency by detecting and quantifying relative bacterial load of P. gingivalis within saliva samples collected from clinical subjects. As a result of gRT-PCR, we confirmed that relative bacterial load of *P. gingivalis* was detected and guantified within all samples of positive control and periodontitis groups. On the contrary, negative results were confirmed in both negative control and healthy groups. Additionally, as a result of comparison with next-generation sequencing (NGS)-based 16S metagenome profiling data, we confirmed relative bacterial load of P. gingivalis, which was not identified on bacterial classification table created through 16S microbiome analysis, in gRT-PCR results. It showed that an approach to quantifying specific microorganisms by applying qRT-PCR method could solve microbial misclassification issues at species level of an NGS-based 16S microbiome study. In this respect, we suggest that *P. gingivalis*-specific primer/probe set introduced in present study has efficient applicability in various oral healthcare industries, including periodontitis-related microbial molecular diagnosis field.

Keywords: molecular diagnosis, *Porphyromonas gingivalis*, quantitative real-time PCR, *rgpA* gene

Availability: Microbiome data (relative abundance table) is exported in QIIME2 microbiome analysis pipeline using python 3.9.7, and is available at https://qiime2.org/.

Introduction

Periodontal disease is one of the most common oral diseases that occur with inflammation inside or outside the tissues surrounding the dental structure (e.g., gums, gingival, and periodontal site), and it has the potential to cause other oral diseases, such as tooth loss [1,2]. In particular, *Porphyromonas gingivalis* is a gram-negative anaerobic microorganism that inhabits the subgingival plaque or saliva, and is classified as a component of periodontitis-associated core pathogenic bacteria (red complex) that affects host immune response and causes chronic periodontitis disease. Therefore, many oral healthcare-related industries have focused on an accurate diagnosis of *P. gingivalis* and research to inhibit intraoral colonization [3].

The quantitative real-time polymerase chain reaction (qRT-PCR) method is generally used in the direct to customers-based molecular diagnostic fields as a powerful tool for quickly and accurately diagnosing infected or not of certain pathogenic microorganisms [4,5]. Although the next-generation sequencing (NGS)based metagenome sequencing method (e.g., 16S microbiome or shotgun metagenome sequencing methods) is possible to determine various microbial communities within each body site associated with human diseases, the qRT-PCR method is still preferred for researchers to detect and quantify specific microorganisms due to the time-cost burden of experiments and analysis of NGS technology [6,7]. Moreover, the high fluorescence sensitivity of qRT-PCR equipment is optimized for specific gene detection and quantification, making it more efficient to manufacture and commercialize a qRT-PCR-based microbial diagnostic kit that detects particular pathogens [8,9]. Although various molecular diagnostic kits or analytic platforms have been launched based on qRT-PCR technology, the market for microbial molecular diagnosis related to periodontal disease is still small compared to other similar fields such as the probiotic industry based on the gut microbiome study. One of the various reasons related to it is that the clinical sampling method performed to diagnose periodontal disease can cause discomfort to consumers in an invasive way that generally requires scratching microbial plaque within the subgingival site and is so difficult to detect in a healthy oral environment [10-12]. However, several recently conducted NGS-based oral microbiome studies have shown that microbial composition detected from saliva or mouthwash can somewhat represent the microbial community in the subgingival plaque sample [10,12]. In particular, Kim et al. [12] showed through the comprehensive 16S microbiome study that there is a significant difference in the relative bacterial frequency of P. gingivalis within saliva samples between periodontitis patients and healthy subjects.

In this study, we designed a specific primer and probe sets for detecting P. gingivalis within saliva samples collected from periodontitis patient (four adult per each group) and healthy groups through qRT-PCR method. The specific primer and probe sets were designed to target the coding region of the bacterial virulence factor 'rgpA' gene (or arg-gingipain A), which has been founded to have a high-expression rate during the biofilm formation process of *P. gingivalis* [13]. We showed that the difference in relative bacterial load of P. gingivalis between each comparison group, which was detected and quantified through the qRT-PCR method, is similar to the NGS-based 16S V3-V4 metagenome sequencing results performed using the same clinical sample. Moreover, we cross-validated the species-specificity of the primer and probe sets through comparative experiments using positive control (two strains belonging to P. gingivalis species) and negative control (genomic DNA pools including other bacterial strains). Therefore, through this application report, we demonstrate that the P. gingiva*lis*-specific primer and probe sets designed in this study have been confirmed high species-specificity that may be applied to the microbial molecular diagnosis of periodontal disease.

Methods and Results

P. gingivalis-specific primer and probe sets design method

First, we designed a species-specific primer and probe sets for detecting and quantifying the bacterial load of *P. gingivalis* within human saliva samples through qRT-PCR, and the process was as follows.

Target gene selection

The *rgpA* is a *P. gingivalis*–specific gene, one of the various virulence gene complexes of *P. gingivalis*, known to increase the risk of periodontal disease by inhibiting the host immune defense system and contributing to the initial inflammatory response and host oral tissue damage around the gums [14].

Obtain of coding sequence information

We obtained the sequence information (FASTA format) of the *rgpA* gene coding sequence region (CDS) for 22 different *P. gingivalis* strains (at the strain level) on the National Center for Biotechnology Information (NCBI DB) reference database (Supplementary Table 1).

Selection of target-specific sequence regions based multiple sequence alignment method

For the selection of primer sequence that can detect all 22 different P.

	A_F1 4530	4540	4550	rgpA_I		4590	rgpA_I	R1 4
1 NC 01072	· · · · · · · · · · · · · · · · · · ·	Variab	ole	///////////////////////////////////////		Consist		GAGANGGA
I NZ CP011 I NZ CP012 I NZ CP024 Reg	gion 🛄	Regio	n			Regio	n ::::	
1 NZ CP024								
1 NZ CP024 1 NZ CP024	GG.	.CTTA	CCA. A					
1 NZ CP024				A 				•••••
1 NZ CP024								
1 NZ CP024	G	.CTTA	CCAA					
1 NZ CP025								
1 N2 CP025 1 N2 CP073 1 N2 CP073		.CTT.	CCAA					
1 N2_CP073	gnments Taxonomy	/						<u></u>
Descriptions Graphic Summary Alig	gnments Taxonomy							<u></u>
Descriptions Graphic Summary Alig Reports Lineage Organism Tax		Score	Number of Hits		Description			<u></u>
Descriptions Graphic Summary Alig Reports Lineage Crganism Tax 12 sequences selecter ?	konomy		Number of Hits	C	Description			<u></u>
Descriptions Graphic Summary Alig Reports Lineage Crganism Tax 42 sequences selecter Organism Porphyromonas . Porphyromonas gingivalis	Blast Name CFB.group bacteria CFB.group bacteria		<u>42</u> <u>35</u>		•			
Descriptions Graphic Summary Alig Reports Lineage Crganism Tax 42 sequences selecter O Organism Porphyromonas • Porphyromonas gingivalis • Porphyromonas gingivalis ATCC 33277	Blast Name CFB group bacteria CFB group bacteria CFB group bacteria	Score 196	42 35 2	Porphyromor	nas gingivalis ATC			
Descriptions Graphic Summary Alig Reports Lineage Crganism Tax 42 sequences selected Organism Porphyromonas gingivalis - Porphyromonas gingivalis ATCC 33277 - Porphyromonas gingivalis ATA1-28	Blast Name CFB. group bacteria CFB. group bacteria CFB. group bacteria CFB. group bacteria	Score 196 196	<u>42</u> <u>35</u>	Porphyromor Porphyromor	nas gingivalis ATO	A1-28 hit		
Descriptions Graphic Summary Alig Reports Lineage Organism Tax 12 sequences selected O Organism Porphyromonas - Porphyromonas gingivalis - Porphyromonas gingivalis - Porphyromonas gingivalis ATA1-28 - Porphyromonas gingivalis 381	Blast Name CFB group bacteria CFB group bacteria CFB group bacteria CFB group bacteria CFB group bacteria	Score 196 196 196	42 35 2	Porphyromor Porphyromor Porphyromor	nas gingivalis AT(nas gingivalis A7/ nas gingivalis 381	<u>A1-28 hit</u> <u>1 hits</u>		
Descriptions Graphic Summary Alig teports Lineage Crganism Tax 2 sequences selected Organism Porphyromonas - Porphyromonas, gingivalis - Porphyromonas, gingivalis - Porphyromonas, gingivalis ATA1:28 - Porphyromonas, gingivalis 381 - Porphyromonas, gingivalis 381 - Porphyromonas, gingivalis AJW4	Blast Name CFB. group bacteria CFB. group bacteria CFB. group bacteria CFB. group bacteria CFB. group bacteria CFB. group bacteria	Score 196 196	42 35 2	Porphyromor Porphyromor Porphyromor Porphyromor	nas gingivalis ATC nas gingivalis A7/ nas gingivalis 387 nas gingivalis AJN	<u>A1-28 hit</u> <u>1 hits</u> <u>N4 hits</u>		
Descriptions Graphic Summary Alig teports Lineage Organism Tax 2 dequences selected O Organism Parphyromonas gingivalis - Porphyromonas gingivalis - Porphyromonas gingivalis ATAL-28 - Porphyromonas gingivalis 381	Blast Name CFB group bacteria CFB group bacteria CFB group bacteria CFB group bacteria CFB group bacteria	Score 196 196 196 196	42 35 2 1 1 1	Porphyromor Porphyromor Porphyromor Porphyromor Porphyromor	nas gingivalis AT(nas gingivalis A7/ nas gingivalis 381	A1-28 hit <u>1 hits</u> <u>N4 hits</u> C60 hits		

Fig. 1. Primer/probe set design process for detection and quantification of relative bacterial load for *Porphyromonas gingivalis*. These figures show the design process of the primer/probe set for quantitative real-time polymerase chain reaction to detect and quantify the relative bacterial load of *Porphyromonas gingivalis*. (A) Multiple sequence alignment results using the Bioedit 7.2 ver. bioinformatics tool show that the primer/probe was designed within consistent regions (highlighted by green box) on the *rgpA* gene of all *P. gingivalis* strains (strain level) annotated on the NCBI database (rgpA_F1: forward primer sequence region; rgpA_R1: reverse primer sequence region; rgpA_P1: probe sequence region). (B) Basic Local Alignment Search Tool (BLAST) test results showing the species-specificity for *P. gingivalis* of primer set.

 Table 1. Sequence information of Porphyromonas gingivalis-specific

 primer and probe sets

Primer/Probe name	Sequence (5' to 3')	Length (mer)
rgpA_Forward	TATCGTTACAGGACAGGGTG	20
rgpA_Reverse	TCCATCTCCTGCGATCCACA	20
rgpA_Probe	56-FAM/TACGAACCCGGAACCTG- CATCC/3SFCQ1	22

gingivalis strains, we performed the multiple sequence alignment (MSA) of each CDS information using the BioEdit 7.2v software (Fig. 1A). We selected two consistent regions confirmed through the MSA method as forward and reverse primer sequences (Table 1).

In silico test

We pre-validated the suitability of the experimental application (Tm value, GC%, and potential for primer-dimer to form, etc.) about the selected primer pair using the 'Oligo calc (http://biotools.nubic.northwetern.edu/)' and 'Oligo Analysis (http://

www.operon.com/tools/oligo-analtsis-tool.aspx)' open web tool. Next, we confirmed the *P. gingivalis*-specificity about each primer sequence and amplification region contained in each primer pair through the NCBI Nucleotide BLAST (https://blast.ncbi.nlm.nih. gov) (Fig. 1B).

Specific-probe design

Finally, we designed a specific probe in the region between each selected primer sequence (Table 1, Fig. 1A). The fluorescent reporter dye and quencher applied to the probe design are 6-FAM (6-carboxyfluorescein) and SFCQ1 (SFC probe, South Korea), which were attached to the *S*' and 3' end regions on the selected probe sequence, respectively.

Comprehensive validation for primer/probe-specificity using qRT-PCR method

We set up the comparison groups to verify experimental detection accuracy of the pre-designed *P. gingivalis*-specific primer/probe

set (Table 2). A total of eight people (four periodontitis and four healthy oral groups) participated in this study, and we obtained saliva samples from them for identifying the bacterial load of *P. gingivalis* within oral cavity. All clinical experiments (recruitment of participants and collection of clinical samples) conducted for this study were approved by the Institutional Review Board of Dankook University Hospital (IRB numbers: 2020-10-015). All clinical examinations were performed by a dentist, who measured pocket depth, clinical attachment loss, gingival index, and plaque index through the full arch (Table 3). All collected clinical specimens were used to extract metagenomic DNA (mDNA) for polymerase chain reaction (PCR), qRT-PCR, and 16S microbiome analysis. The mDNA was extracted using the QIAamp PowerFecal Pro DNA Kit (QIAGEN, Hilden, Germany), and all experimental processes were performed according to the formal protocol guide pro-

 Table 3. Demographic characteristic and clinical parameters of periodontal disease group and healthy group

Characteristic	Healthy	Periodontitis	p-value
Sex			
Male	3	4	
Female	1	0	
Age (y) (mean ± SD)	28.1 ± 10.55	38.8 ± 13.40	< 0.01
Examination			
PD (mm)	2.49 ± 0.51	5.97 ± 2.27	< 0.01
CAL (mm)	2.31 ± 0.22	6.12 ± 2.09	< 0.01
GI	0.17 ± 0.46	1.69 ± 0.94	< 0.01
PI	0.35 ± 0.62	1.28 ± 0.92	< 0.01

SD, standard deviation; PD, probing depth; CAL, clinical attachment loss; GI, gingival index; PI, plaque index

Table 2. Sample information applied to this study

vided in the Kit. The 'ZymoBIOMICS standard DNA', which belongs to the negative control group, is a mixed DNA sample in which the genomic DNAs of 10 different bacterial species (*Listeria* monocytogenes, Pseudomonas aeruginosa, Bacillus subtilis, Saccharomyces cerevisiae, Escherichia coli, Salmonella enterica, Lactobacillus fermentum, Enterococcus faecalis, Cryptococcus neoformans, and Staphylococcus aureus) are at a uniform rate (Zymo Research, Irvine, CA, USA). The experimental validation process for primer/ probe activity is as follows.

PCR validation for confirming primer specificity

First, we performed PCR validation to confirm whether the pre-designed primer binds specifically to the *rgpA* gene CDS region on the *P. gingivalis* genome (Fig. 2). As a result of PCR validation, we confirmed that a correct target amplicon DNA band (approximately 110 bp) was located within samples corresponding to positive control and periodontitis groups. In contrast, amplicon DNA was not identified within samples corresponding to the negative control and healthy group.

Comparison of bacterial load about P. gingivalis within each sample using qRT-PCR

Next, we performed qRT-PCR to confirm the relative bacterial load of the *P. gingivalis* within each comparison group (Table 4, Fig. 3, Supplementary Table 2). Prior to the experiment, the concentration of all double-strand DNA samples was checked through the Qubit Fluorometer 4.0 v at the Center for Bio-medical Engineering Core Facility (Dankook University, South Korea) and 1× dsDNA HS Assay kit (Thermo Fisher Scientific, Waltham, MA, USA), and then normalized to the consistent concentration (10

Groups	Sample name	DNA type	Extracted raw DNA concentration (ng/ μ L)
Positive control	Porphyromonas gingivalis ATCC 49417	Single species DNA	13.50
	P. gingivalis ATCC 33277		4.00
Periodontitis	Case_01	Metagenomic DNA	46.00
	Case_02		47.20
	Case_03		26.60
	Case_04		15.30
Healthy	Normal_01	Metagenomic DNA	47.20
	Normal_02		33.20
	Normal_03		29.80
	Normal_04		21.20
Negative control	Epidermidibacterium keratini EPI-7 [™]	Single species DNA	10.90
	ZymoBIOMICS standard DNA	Multiple species DNA	10.00
	NTC (DW)	None DNA	-




Fig. 2. Polymerase chain reaction (PCR) validation result for evaluating primer binding activity of *Porphyromonas gingivalis*-specific primer set. Gel-electrophoresis showing the results of evaluating primer binding activity for the rgpA gene region on the *P. gingivalis* genome through PCR validation.

		Ct value of each dilutied DNA sample							
Groups	Samples	10 ng (raw co	oncentration)	1 ng (10 ⁻¹ dilution)		0.1 ng (10 ⁻² dilution)		0.01 ng (10^{-3} dilution)	
		Avg	SD	Avg	SD	Avg	SD	Avg	SD
Positive control	P. gingivalis ATCC 49417	19.67	0.04	22.34	0.18	25.12	0.04	27.85	0.04
	P. gingivalis ATCC 33277	19.45	0.06	22.02	0.25	24.8	0.01	27.33	0.02
Periodontitis	Case_01	27.28	0.09	30.33	0.06	33.45	0.11	37.04	0.01
	Case_02	27.49	0.02	30.18	0.19	33.13	0.03	35.51	0.12
	Case_03	26.65	0.04	30.02	0.15	32.57	0.10	35.01	0.19
	Case_04	31.01	0.18	33.31	0.05	36.98	0.79	39.21	0.08
Healthy	Normal_01	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Normal_02	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Normal_03	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Normal_04	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Negative control	Epidermidibacterium keratini EPI-7T	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	ZymoBIOMICS standards	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	NTC (D.W.)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

Table 4. Comparison of bacterial load about *Porphyromonas gingivalis* using gRT-PCR

qRT-PCR, quantitative real-time polymerase chain reaction; Avg, average Ct value for each triplicate calculations; SD, standard deviation value for each triplicate calculations; N/A, no detection.

 $ng/\mu L$). This process is essential to compare the Ct value (cycle threshold) reflecting the relative frequency of *P. gingivalis* on the same mDNA concentration within saliva samples of each clinical subject (healthy and periodontitis). The same experimental conditions were set by normalizing the other microbial genomic DNA concentrations (positive and negative controls) equally with the clinical comparison group. After normalizing all DNA concentrations, to confirm the relative bacterial load of *P. gingivalis* as various DNA concentrations range, we performed qRT-PCR by diluting each template DNA (10^1 , 10^0 , 10^{-1} , 10^{-2}) through 10-fold serial dilution method (using CFX Opus 96 Real-Time PCR System, Bio-

Rad, Hercules, CA, USA; The running condition of the qRT-PCR is as follows: pre-denaturation 95.0°C, 10 min; denaturation 95.0°C, 10 s; annealing 64.0°C, 15 s; elongation/flour detection 72.0°C, 15 s; and total PCR cycle was 45). As a result of qRT-PCR, we confirmed through the Ct value measurement results that the *P. gingivalis rgpA* gene was detected from all diluted DNA samples within the positive control and periodontitis groups. In contrast, *P. gingivalis* detection through confirming *rgpA* gene amplification could not be confirmed within the healthy and negative control groups (denoted 'N/A' in Table 4 and Supplementary Table 2). Although high-average Ct values (31.01_10 ng; 33.31_1 ng;




Fig. 3. Quantitative real-time polymerase chain reaction (qRT-PCR) result reflecting the relative bacterial load of *Porphyromonas gingivalis* within each DNA sample. The amplification plot shows whether the relative bacterial load of *P. gingivalis* was detected within all samples in each comparison group. The x-axis on the plot indicates the running cycle number of qRT-PCR, and the y-axis indicates relative fluorescence units (RFU).

36.98_0.1 ng; 39.21_0.01 ng) were calculated from one sample in the periodontitis group (Case_04), our data show that *P. gingivalis* with low relative abundance within saliva sample can also be detected via qRT-PCR.

Comparison of bacterial detection between qRT-PCR and NGSmetagenome profiling

Finally, we checked the relative bacterial frequency of P. gingivalis within all samples of the healthy and periodontitis groups (four samples per each group) through the NGS-based 16S metagenome profiling analysis using the QIIEM2 microbiome analysis pipeline and compared it with the qRT-PCR result (Supplementary Table 3). The 16S microbiome analysis about eight clinical samples was performed through bacterial 16S V3-V4 amplicon sequencing using pre-extracted mDNA samples (Table 4), and the Illumina Miseq sequencer platform (300PE sequencing, Illumina, San Diego, CA, USA) was applied to the generation of raw sequencing data. We analyzed the paired-end 16S V3-V4 raw sequencing read files (i.e., fastq.gz files) using two bioinformatics pipelines: QIIME2 and DADA2 (version 1.2.0). The executing software for QIIME2 was Python software and for DADA2 was R (version 3.3.2). First, the sequencing reads were sorted using unique barcodes for each PCR product. The barcode, linker, and primer sequences were then removed from the original sequencing reads. And the removed reads were merged paired-end reads using FLASH v 1.2.11. The merged reads containing two or more ambiguous nucleotides, those with a low-quality score (average

score < 20), or reads shorter than 300 bp, were filtered out. All pre-filtered sequence data quality control and feature table construction were completed using the Divisive Amplicon Denoising Algorithm 2 (DADA2) in the QIIME 2 plugin, which detects and corrects amplicon errors and filters out PhiX chimeric sequences. The pre-processed reads from each sample were used to calculate the number of Amplicon Sequence Variants (ASVs). The number of ASVs was determined by clustering the sequences from each sample using a 99% sequence identity cut-off using DADA2 pipeline-based QIIME2 software. Taxonomic relative abundance was counted with Naïve Bayes classifier using a confidence threshold of 0.7% derived from the ASVs data for each sample. The classified microbial composition based on SILVA v138 16S rRNA gene database was normalized using the value calculated from the taxonomy abundance count divided by the number of pre-processed reads for each sample. The relative bacterial frequency of P. gingivalis in each group identified through microbiome analysis was as follows: Healthy group (Normal 01: 0.00, Normal 02: 0.00, Normal 03: 0.00, and Normal 04: 0.00); Periodontitis group (Case 01: 0.00, Case 02: 0.16, Case 03: 0.16, and Case 04: 0.00). As with the qRT-PCR results, we confirmed that *P. gingivalis* was not detected within all samples of the healthy group. However, unlike the qRT-PCR results, we identified that the relative bacterial frequency of *P. gingivalis* was not calculated within two samples (Case 01 and Case 04) of the periodontitis group. These results demonstrated that qRT-PCR also sensitively captured the low abundance for P. gingivalis, which was not confirmed by the NGS-

based metagenome sequencing method.

Discussion

We introduced a high-sensitivity primer/probe set that can specifically detect and quantify the relative bacterial load of periodontitis-related pathogenic bacteria 'Porphyromonas gingivalis' within human saliva sample. This primer/probe set was designed to target the P. gingivalis-specific gene, rgpA gene, and has never been applied as a target gene to diagnose the P. gingivalis. Moreover, we suggest that the microbial detection accuracy of our primer/probe set is excellent because we designed a primer/probe set to cover the rgpA gene of all P. gingivalis strains (strain level) annotated on the NCBI reference database. By comparing bacterial detection and quantification results for P. gingivalis between NGS-based metagenome profiling data and qRT-PCR, we cross-validated that the qRT-PCR method using the *rgpA* gene-specific primer/ probe set showed higher detection sensitivity than the microbiome analysis targeting the bacterial 16S ribosomal RNA gene. This result shows that microbial misclassification problem at the species level due to high sequence similarity within the hyper-variable regions on the bacterial 16S rRNA gene occurring from 16S microbiome analysis can be overcome through the qRT-PCR method [15,16]. In this respect, we suggest that the *P. gingivalis*-specific primer/ probe set introduced in the present study has efficient applicability in various oral healthcare industries, including the periodontitis-related microbial molecular diagnosis field.

ORCID

Jinuk Jeong: https://orcid.org/0000-0002-7283-0949 Yunseok Oh: https://orcid.org/0000-0002-6957-2715 Junhyeon Jeon: https://orcid.org/0000-0002-6549-3243 Dong-Heon Baek: https://orcid.org/0000-0002-9450-4247 Dong Hee Kim: https://orcid.org/0000-0002-5056-7406 Kornsorn Srikulnath: https://orcid.org/0000-0002-5985-7258 Kyudong Han: https://orcid.org/0000-0001-6791-2408

Authors' Contribution

Conceptualization: KH. Data curation: YO, JJ (Junhyeon Jeon), JJ (Jinuk Jeong). Formal analysis: DHK, KS, JJ (Junhyeon Jeon). Methodology: YO, DHK, JJ (Junhyeon Jeon). Writing - original draft: JJ (Junhyeon Jeon), KH. Writing - review & editing: all authors.

Conflicts of Interest

No potential conflict of interest relevant to this article was reported.

Acknowledgments

The authors gratefully acknowledge Center for Bio-Medical Engineering Core Facility at Dankook University for providing critical reagents and equipment. The research institute has been supported by the VIP system as a part of the University Innovation Support Program 2022 of Dankook University. The research institute has been supported by the Department of Microbiology was supported through the Research-Focused Department Promotion & Interdisciplinary Convergence Research Project as a part of the Support Program for University Development for Dankook University in 2022.

Supplementary Materials

Supplementary data can be found with this article online at http://www.genominfo.org.

References

- 1. Abusleme L, Hoare A, Hong BY, Diaz PI. Microbial signatures of health, gingivitis, and periodontitis. Periodontol 2000 2021; 86:57-78.
- Kononen E, Gursoy M, Gursoy UK. Periodontitis: a multifaceted disease of tooth-supporting tissues. J Clin Med 2019;8:1135.
- **3.** Damgaard C, Danielsen AK, Enevold C, Massarenti L, Nielsen CH, Holmstrup P, et al. *Porphyromonas gingivalis* in saliva associates with chronic and aggressive periodontitis. J Oral Microbiol 2019;11:1653123.
- 4. Esteves LM, Bulhoes SM, Branco CC, Carreira T, Vieira ML, Gomes-Solecki M, et al. Diagnosis of human leptospirosis in a clinical setting: real-time PCR high resolution melting analysis for detection of leptospira at the onset of disease. Sci Rep 2018;8: 9213.
- Maurin M, Hammer L, Gestin B, Timsit JF, Rogeaux O, Delavena F, et al. Quantitative real-time PCR tests for diagnostic and prognostic purposes in cases of legionellosis. Clin Microbiol Infect 2010;16:379-384.
- 6. Jeong J, Mun S, Oh Y, Cho CS, Yun K, Ahn Y, et al. A qRT-PCR method capable of quantifying specific microorganisms compared to NGS-based metagenome profiling data. Microorgan

isms 2022;10:324.

- 7. O'Sullivan DJ, Giblin L, McSweeney PL, Sheehan JJ, Cotter PD. Nucleic acid-based approaches to investigate microbial-related cheese quality defects. Front Microbiol 2013;4:1.
- 8. Ambrose JH, Sekaran SD, Azizan A. Dengue virus NS1 protein as a diagnostic marker: commercially available ELISA and comparison to qRT-PCR and serological diagnostic assays currently used by the State of Florida. J Trop Med 2017;2017:8072491.
- 9. Qi Y, Cui L, Ge Y, Shi Z, Zhao K, Guo X, et al. Altered serum microRNAs as biomarkers for the early diagnosis of pulmonary tuberculosis infection. BMC Infect Dis 2012;12:384.
- Boutaga K, Savelkoul PH, Winkel EG, van Winkelhoff AJ. Comparison of subgingival bacterial sampling with oral lavage for detection and quantification of periodontal pathogens by real-time polymerase chain reaction. J Periodontol 2007;78:79-86.
- Boutaga K, van Winkelhoff AJ, Vandenbroucke-Grauls CM, Savelkoul PH. Comparison of real-time PCR and culture for detection of *Porphyromonas gingivalis* in subgingival plaque samples. J Clin Microbiol 2003;41:4950-4954.

- Kim YT, Jeong J, Mun S, Yun K, Han K, Jeong SN. Comparison of the oral microbial composition between healthy individuals and periodontitis patients in different oral sampling sites using 16S metagenome profiling. J Periodontal Implant Sci 2022; 52:394-410.
- 13. Liang G, Shi H, Qi Y, Li J, Jing A, Liu Q, et al. Specific anti-biofilm activity of carbon quantum dots by destroying *P. gingivalis* biofilm related genes. Int J Nanomedicine 2020;15:5473-5489.
- 14. Zhang R, Yang J, Wu J, Sun WB, Liu Y. Effect of deletion of the rgpA gene on selected virulence of *Porphyromonas gingivalis*. J Dent Sci 2016;11:279-286.
- Callahan BJ, Grinevich D, Thakur S, Balamotis MA, Yehezkel TB. Ultra-accurate microbial amplicon sequencing with synthetic long reads. Microbiome 2021;9:130.
- 16. Jeong J, Yun K, Mun S, Chung WH, Choi SY, Nam YD, et al. The effect of taxonomic classification by full-length 16S rRNA sequencing with a synthetic long-read technology. Sci Rep 2021; 11:1727.